

# ECG 心电(G)

## ECG CABLES&LEADWIRES 心电导联线


- Fixed One-piece, Trunk cable, Leadwire types available  
一体化，主电缆，导联线可用
- Latex-free  
不含乳胶
- 3 leads, 5 leads, Snap, Pinch available  
3 导，5 导，扣式，夹式可用
- AHA and IEC color coding available  
美标，欧标可用
- Different types are Compatible with most of patient monitors  
多款样式的探头兼容大多数病人监护仪
- Easy to clean and repair  
容易清洗和维修




Product types 产品类别	Photo 图片	P/N Example 型号示例	P/N Guide 型号识别
Fixed one-piece ECG cable 一体化心电导联线		<u><a href="#">G 5 1 01 S</a></u>	<b>G:</b> ECG cables <b>5:</b> 5 leads (3: 3 leads) <b>1:</b> AHA (2: IEC) <b>01:</b> General (10: GE-DATEX; 12: GE-Marquette) <b>S:</b> fixed Snap (P: Fixed pinch)
Trunk cable 心电主电缆		<u><a href="#">G 5 2 12 MQ</a></u>	<b>G:</b> ECG cables <b>5:</b> 5 leads (3: 3 leads) <b>2:</b> IEC (1: AHA) <b>12:</b> GE-Marquette (10: GE-DATEX) <b>MQ:</b> MQ type
Leadwire 心电导联线		<u><a href="#">G 3 1 1 AA</a></u>	<b>G:</b> ECG cables <b>3:</b> 3 leads (5: 5 leads) <b>1:</b> AHA (2: IEC) <b>1:</b> Pinch (2: snap) <b>AA:</b> AA type
<b>Warranty:</b> six-month warranty for ECG cables, trunk cables and leadwires <b>Package:</b> 1pcs/bag <b>保修:</b> 6 个月 <b>包装:</b> 1 个/包			

Monitor Manufacturers 监护仪厂家	Leads/label 导联数	Trunk cable 主电缆	Trunk Cable yoke# 分线盒	Leadwire-Pinch 夹式导联线	Leadwire-Snap 扣式导联线	Leadwire Conn.# 导联线接头	Fixed-Pinch 夹式一体化	Fixed-Snap 扣式一体化	ECG Conn. 心电接头
<b>General 6 pins 普通 6 针</b>									
Straight connector	3/AHA	G3140DN	CN151	G311DN	G312DN	CN153	G3140P	G3140S	CN101B
	3/IEC	G3240DN		G321DN	G322DN		G3240P	G3240S	
	5/AHA	G5140DN	CN152	G511DN	G512DN		G5140P	G5140S	
	5/IEC	G5240DN		G521DN	G522DN		G5240P	G5240S	
Angle connector	3/AHA	G3101DN	CN151	G311DN	G312DN	CN153	G3101P	G3101S	CN101
	3/IEC	G3201DN		G321DN	G322DN		G3201P	G3201S	
	5/AHA	G5101DN	CN152	G511DN	G512DN		G5101P	G5101S	
	5/IEC	G5201DN		G521DN	G522DN		G5201P	G5201S	
LL style	3/AHA	G3140LL	CN154	G311LL	G312LL	CN156	G3101P	G3101S	CN101B
	3/IEC	G3240LL		G321LL	G322LL		G3201P	G3201S	
	5/AHA	G5140LL	CN155	G511LL	G512LL		G5101P	G5101S	
	5/IEC	G5240LL		G521LL	G522LL		G5201P	G5201S	
<b>General 6 pins with Resistor 普通 6 针带电阻</b>									
Normal 6 pins conn. With 1k resistor, Straight connector	3/AHA	G3141DN	CN151	G311DN	G312DN	CN153	G3141P	G3141S	CN101B
	3/IEC	G3241DN		G321DN	G322DN		G3241P	G3241S	
	5/AHA	G5141DN	CN152	G511DN	G512DN		G5141P	G5141S	
	5/IEC	G5241DN		G521DN	G522DN		G5241P	G5241S	
<b>Artema-S/W 阿提玛</b>									
Series 8000, Athena,Diascope 1-2	3/AHA						G3103P	G3103S	CN102
	3/IEC						G3203P	G3203S	
<b>Bionet</b>									
BM3	3/AHA						G3104P	G3104S	CN103
	3/IEC						G3204P	G3204S	
	5/AHA						G5104P	G5104S	
	5/IEC						G5204P	G5204S	
<b>Biosys</b>									
	3/AHA						G3105P	G3105S	CP6
	3/IEC						G3205P	G3205S	
	5/AHA						G5105P	G5105S	
	5/IEC						G5205P	G5205S	
<b>Colin 科林</b>									
BP 88 - 306 - 508 Plastic connector	3/AHA						G3106P	G3106S	CN104
	3/IEC						G3206P	G3206S	
	5/AHA						G5106P	G5106S	
	5/IEC						G5206P	G5206S	



Monitor Manufacturers 监护仪厂家	Leads/label 导联数	Trunk cable 主电缆	Trunk Cable yoke# 分线盒	Leadwire-Pinch 夹式导联线	Leadwire-Snap 扣式导联线	Leadwire Conn.# 导联线接头	Fixed-Pinch 夹式一体化	Fixed-Snap 扣式一体化	ECG Conn. 心电接头
<b>Colin 科林</b>									
Mental L connector 	3/AHA						G3107P	G3107S	CN105
	3/IEC						G3207P	G3207S	
	5/AHA						G5107P	G5107S	
	5/IEC						G5207P	G5207S	
<b>Drager-Siemens 德尔格-西门子</b>									
Sirecust 341-342 R Series 200-400-600-700- 900-1200	3/AHA	G3108DR	CN159	G311DR	G312DR	CN161	G3108P	G3108S	CN106
	3/IEC	G3208DR		G321DR	G322DR		G3208P	G3208S	
	5/AHA	G5108DR	CN160	G511DR	G512DR	CN162	G5108P	G5108S	
				5956490					
	5/IEC	G5208DR		G521DR	G522DR		G5208P	G5208S	
			5956508						
SC5000, SC6000, SC7000, SC8000, SC9000 	3/AHA						G3131P	G3131S	CN128
	3/IEC						G3231P	G3231S	
<b>Fukuda Denshi 福田</b>									
DS-7000 DS-7100 DS-7200 DS-7300	5/AHA	G5109DR	CN160	G511DR	G512DR	CN162	G5109P	G5109S	CN107
	5/IEC	G5209DR		G521DR	G522DR		G5209P	G5209S	
DS-5300W DS-5100E DS-7001 	3/AHA						G3133P	G3133S	CN123
	3/IEC						G3233P	G3233S	
	5/AHA						G5133P	G5133S	
	5/IEC						G5233P	G5233S	
<b>GE-Critikon 强生</b>									
Series 8700 - Dinamap + Dinamap MPS Select Dinamap MPS Select Portable	3/AHA	G3102DX	CN163	G311DX	G312DX	CN165	G3102P	G3102S	CN101
	3/IEC	G3202DX		G321DX	G322DX		G3202P	G3202S	
	5/AHA	G5102DX	CN164	G511DX	G512DX	CN166	G5102P	G5102S	
	5/IEC	G5202DX		G521DX	G522DX		G5202P	G5202S	
<b>GE-Datex Ohmeda 德恩</b>									
AS/3, CS/3, Cardiocap 1, Cardiocap 2, Cardiocap 5, S/5, Satellite Plus	3/AHA	G3110DX	CN163	G311DX	G312DX	CN165	G3110P	G3110S	CN125
		545302		8001958	545327				
		545307		545317					
	3/IEC	G3210DX		G321DX	G322DX		G3210P	G3210S	
		545305 545300	8001960						
		545304	545315						
	5/AHA	G5110DX	CN164	G511DX	G512DX	CN166	G5110P	G5110S	
		545303		8001959	545328				
		545308		545318					
	5/IEC	G5210DX		G521DX	G522DX		G5210P	G5210S	
	545301	8001961							
	545306	545316							




Monitor Manufacturers 监护仪厂家	Leads/label 导联数	Trunk cable 主电缆	Trunk Cable yoke# 分线盒	Leadwire-Pinch 夹式导联线	Leadwire-Snap 扣式导联线	Leadwire Conn.# 导联线接头	Fixed-Pinch 夹式一体化	Fixed-Snap 扣式一体化	ECG Conn. 心电接头
<b>GE-Hellige 海力格</b>									
Eagle 1000, Def. SCP 900,Cardioserv, SMC,SMS,SMK,SMU,SMV	3/AHA	G3111DX	CN163	G311DX	G312DX	CN165	G3111P	G3111S	CN108
	3/IEC	G3211DX		G321DX	G322DX		G3211P	G3211S	
	5/AHA	G5111DX	CN164	G511DX	G512DX	CN166	G5111P	G5111S	
	5/IEC	G5211DX		G521DX	G522DX		G5211P	G5211S	
<b>GE-Marquette 马奎</b>									
Dash 2000/3000/4000, Dinamap PRO 1000,Eagle Systems, MAC-Lab,Solar 3000/8000/ 9500,Tram 100/200/300/400/ 500/600/800/Tram AR,Unity Systems	3/AHA	G5112MQ	CN167	G311MQ	G312MQ	CN168	G3112P	G3112S	CN109
		412931-001		412682-001	411203-001				
		412931-021		412682-002	411203-002				
	3/IEC	G5212MQ		G321MQ	G322MQ		G3212P	G3212S	
		2017003-003		412682-003	411203-003				
		412931-002		412682-004	411203-004				
	5/AHA	G5112MQ		G511MQ	G512MQ		G5112P	G5112S	
				412681-001	412681-003				
				412681-002	412681-004				
	5/IEC	G5212MQ		G521MQ	G522MQ		G5212P	G5212S	
10/AHA	G1112MQ								
10/IEC	G1212MQ								
<b>Huntleigh Healthcare</b>									
Ⓝ							G3142P	G5142P	CP8
							G3242P	G5242P	
<b>Kontron 康强</b>									
Micromon, Minimon,Vitalmon, Supermon, Kolormon	3/AHA						G3113P	G3113S	CN110
	3/IEC						G3213P	G3213S	
105, 107, 108, 125, 128, VSM510 Old generation	3/AHA						G3114P	G3114S	CN36
	3/IEC						G3214P	G3214S	
<b>Medtronic-Physio Control 美敦力</b>									
Lifepak5,6,7,8,9(See General(Resistor))									
Lifepak 15, 20 Ⓝ	3/AHA						G3115P	G3115S	CN111
	3/IEC						G3215P	G3215S	
	5/AHA						G5116P	G5115S	
	5/IEC						G5215P	G5215S	
	10/AHA						G1115P	G1115S	CN111
	10/IEC						G1215P	G1215S	



Monitor Manufacturers 监护仪厂家	Leads/label 导联数	Trunk cable 主电缆	Trunk Cable yoke# 分线盒	Leadwire-Pinch 夹式导联线	Leadwire-Snap 扣式导联线	Leadwire Conn.# 导联线接头	Fixed-Pinch 夹式一体化	Fixed-Snap 扣式一体化	ECG Conn. 心电接头
<b>Mennen 美能</b>									
Series 9000-1100-2000	3/AHA	G3116MN	CN173	G311MN	G312MN	CN175	G3116P	G3116S	CN112
	3/IEC	G3216MN		G321MN	G322MN		G3216P	G3216S	
	5/AHA	G5116MN	CN174	G511MN	G512MN	CN176	G5116P	G5116S	
	5/IEC	G5216MN		G521MN	G522MN		G5216P	G5216S	
Envoy	3/AHA	G3117MN	CN173	G311MN	G312MN	CN175	G3117P	G3117S	CN113
	3/IEC	G3217MN		G321MN	G322MN		G3217P	G3217S	
	5/AHA	G5117MN	CN174	G511MN	G512MN	CN176	G5117P	G5117S	
	5/IEC	G5217MN		G521MN	G522MN		G5217P	G5217S	
<b>Mindray 迈瑞</b>									
6 pins with resistor 带电阻 PM-7000, PM-8000series, PM-9000series, new PM-6000	3/AHA	G3141MD	CN159	G311MD	G312MD	CN161	G3141P	G3141S	CN101B
		0010-30-12257 0010-30-12246		G321MD	G322MD		G3241P	G3241S	
	3/IEC	G3241MD		G321MD	G322MD		G3241P	G3241S	
		0010-30-12247	G311MD	G312MD	G3141P	G3141S			
	5/AHA	G5141MD	CN160	G511MD	G512MD	CN162	G5141P	G5141S	
		0010-30-12244		G521MD	G522MD		G5241P	G5241S	
5/IEC	G5241MD	G521MD		G522MD	G5241P		G5241S		
	0010-30-12245	G311MD	G312MD	G3141P	G3141S				
12 pins, no resistor 12 针不带电阻 PM-5000	3/AHA	G3118MD	CN159	G311MD	G312MD	CN161	G3118P	G3118S	CN117
		0010-30-12380 0010-30-12250		G321MD	G322MD		G3218P	G3218S	
	3/IEC	G3218MD	G321MD	G322MD	G3218P	G3218S			
		0010-30-12251							
12 pins with resistor 12 针带电阻 PM-6000	3/AHA	G3124MD	CN159	G311MD	G312MD	CN161	G3124P	G3124S	CN117
		0010-30-12379		G321MD	G322MD		G3224P	G3224S	
	3/IEC	G3224MD		G321MD	G322MD		G3224P	G3224S	
		0010-30-12379	G311MD	G312MD	G3124P	G3124S			
	5/AHA	G5124MD	CN160	G511MD	G512MD	CN162	G5124P	G5124S	
		0010-30-12260		G521MD	G522MD		G5224P	G5224S	
5/IEC	G5224MD	G521MD		G522MD	G5224P		G5224S		
	0010-30-12260								



Monitor Manufacturers 监护仪厂家	Leads/label 导联数	Trunk cable 主电缆	Trunk Cable yoke# 分线盒	Leadwire-Pinch 夹式导联线	Leadwire-Snap 扣式导联线	Leadwire Conn.# 导联线接头	Fixed-Pinch 夹式一体化	Fixed-Snap 扣式一体化	ECG Conn. 心电接头		
<b>Mindray 迈瑞</b>											
12 pins with resistor 12 针带电阻 BENEVIEW T5/T8	3/AHA	G3143MD	CN159	G311MD	G312MD	CN161	G3143P	G3143S	CN127		
		0010-30-42719			0010-30-42726		0010-30-42734				
	3/IEC	G3243MD			G321MD		G322MD			G3243P	G3243S
		0010-30-42719			0010-30-42725		0010-30-42733				
	5/AHA	G5143MD	CN160	G511MD	G512MD	CN162	G5143P	G5143S			
		0010-30-42719			0010-30-42727		0010-30-42735				
5/IEC	G5243MD			G521MD	G522MD			G5243P	G5243S		
	0010-30-42719		0010-30-42728	0010-30-42736							
6 pins , no resistor 不带电阻 MEC1000, MEC1200	3/AHA	G3140MD	CN159	G311MD	G312MD	CN161	G3140P	G3140S	CN101B		
		0010-30-12378 0010-30-12242									
	3/IEC	G3240MD			G321MD		G322MD			G3240P	G3240S
		0010-30-12243									
<b>Mediana</b>											
	3/AHA	G3129DN					G3129P	G3129S	CP8		
	3/IEC	G3229DN					G3229P	G3229S			
<b>MEK</b>											
MP500,MP600,MP1000	3/AHA						G3119P	G3119S	CN30		
	3/IEC						G3219P	G3219S			
	5/AHA						G5119P	G5119S			
	5/IEC						G5219P	G5219S			
6 pins	3/AHA						G3120P	G3120S	CN101B		
	3/IEC						G3220P	G3220S			
	5/AHA						G5120P	G5120S			
	5/IEC						G5220P	G5220S			
<b>Nihon Kohden 光电</b>											
Def. Cardiolife Lifescope serioes 6-8-9-11 6102-6205-6120	3/AHA	G5121NH	CN169	G311NH	G312NH	CN170	G3121P	G3121S	CN114		
	3/IEC	G5221NH			G321NH		G322NH			G3221P	G3221S
	5/AHA	G5121NH			G511NH		G512NH			G5121P	G5121S
	5/IEC	G5221NH			G521NH		G522NH			G5221P	G5221S
BSM-2300 liflescope I,BSM-4100 liflescope P, BSM-5100 liflescope A,BSM-9510liflescope M, BSM-9800 liflescope S	3/AHA	G6122NK	CN171	G311NK	G312NK	CN172	G3122P	G3122S	CN115		
	3/IEC	G6222NK			G321NK		G322NK			G3222P	G3222S
	6/AHA	G6122NK			G611NK		G612NK				
	6/IEC	G6222NK			G621NK		G622NK				
TEC-5521	3/AHA	G3130NH	CN169	G311NH	G312NH	CN170	G3130P	G3130S	CN119		
	3/IEC	G3230NH		G321NH	G322NH		G3230P	G3230S			



Monitor Manufacturers 监护仪厂家	Leads/label 导联数	Trunk cable 主电缆	Trunk Cable yoke# 分线盒	Leadwire-Pinch 夹式导联线	Leadwire-Snap 扣式导联线	Leadwire Conn.# 导联线接头	Fixed-Pinch 夹式一体化	Fixed-Snap 扣式一体化	ECG Conn. 心电接头
<b>Nihon Kohden 光电</b>									
OPV-1500 	3/AHA	G3134NK	CN171	G311NK	G312NK	CN172	G3134P	G3134S	CN126
	3/IEC	G3234NK		G321NK	G322NK		G3234P	G3234S	
<b>Philips-HP 飞利浦</b>									
43100A,43110A,43110M/MC,43120A,43200A, 77000AC,77020A/AC/AR,77025A,77030A,77035A, 77040A,78208A,78210A,78211A,78330A,78331A,78332A ,78333A,78341A,78342A,78345A,78346A,78351T,78352 A,78353A,78353B,78354A,78532A,78532B,78534A/B/C, 78619A,78620A,78660A/B,78661A/B,78670A,78671A,M1 722A/B,M1723A/B,M2400A,M2406A,M2745B,M3500B, M5500B	3/AHA	G3123AA	CN173	G311AA	G312AA	CN175	G3123P	G3123S	CN116
		M1580A		M1603A	M1605A				
	3/IEC	G3223AA	G321AA	G322AA	G3223P	G3223S			
		M1590A	M1613A	M1615A					
78352C,78354C,78834C,862474 (C3), 862475 (C3),862478 (C3),862479 (C3), MP20,MP30,MP40,MP50,MP60,MP70,MP90, M1001A,M1001B,M1002A,M1002B,M1165A, M1166A,M1167A,M1175A,M1177A,M1205A, M1722A,M1722B,M2601A,M2745A,M3000A, M3001A,M3046A,M3500B,M3535A,M3536A, M4375A,M5500B	5/AHA	G5123AA	CN174	G511AA	G512AA	CN176	G5123P	G5123S	CN117
		M1600A		M1623A	M1625A				
5/IEC	G5223AA	G521AA	G522AA	G5223P	G5223S				
	M1610A	M1633A	M1635A						
78352C, 78354C, 78834C / 862472 to 862475 C3 Monitor / 862478, 862479 / 863064, 863065, 863066, 863068 SureSigns VM6, VM8 / M1001A/B, M1002A/B / M3000A, M3001A, M3002A MMS X2 M8102A MP2 Monitor, M8105A MP5 Monitor / M4735A, M3535A/M3536A 	3/AHA	G3124AA	CN173	G311AA	G312AA	CN175	G3124P	G3124S	CN117
		M1500A		M1603A	M1605A				
	3/IEC	G3224AA	G321AA	G322AA	G3224P	G3224S			
		M1510A	M1613A	M1615A					
78352C, 78354C, 78834C / 862472 to 862475 C3 Monitor / 862478, 862479 / 863064, 863065, 863066, 863068 SureSigns VM6, VM8 / M1001A/B, M1002A/B / M3000A, M3001A, M3002A MMS X2 M8102A MP2 Monitor, M8105A MP5 Monitor / M4735A, M3535A/M3536A 	5/AHA	G5124AA	CN174	G511AA	G512AA	CN176	G5124P	G5124S	CN117
		M1520A		M1623A	M1625A				
	5/IEC	G5224AA	G521AA	G522AA	G5224P	G5224S			
		M1530A	M1633A	M1635A					
78352C, 78354C, 78834C / 862472 to 862475 C3 Monitor / 862478, 862479 / 863064, 863065, 863066, 863068 SureSigns VM6, VM8 / M1001A/B, M1002A/B / M3000A, M3001A, M3002A MMS X2 M8102A MP2 Monitor, M8105A MP5 Monitor / M4735A, M3535A/M3536A 	3/IEC	G3224PH	CN179	G311PH	G312PH	CN181			CN117
		M1669A		M1671A	M1673A				
	3/AHA		G321PH	G322PH					
			M1672A	M1674A					
M8102A MP2 Monitor, M8105A MP5 Monitor / M4735A, M3535A/M3536A 	5/AHA	G5224PH	CN180	G511PH	G512PH	CN182			CN117
		M1668A		M1968A	M1644A				
	5/IEC		G521PH	G522PH					
			M1971A	M1645A					
<b>Schiller/Brucker/ODM</b>									
	3/AHA						G3125P	G3125S	CN110
	3/IEC						G3225P	G3225S	
	5/AHA						G5125P	G5125S	
	5/IEC						G5225P	G5225S	
<b>Spacelabs 太空</b>									
Ultraview New generation	3/AHA	G5126SL	CN177	G311SL	G312SL	CN178	G3126P	G3126S	CN118
	3/IEC	G5226SL		G321SL	G322SL		G3226P	G3226S	
	5/AHA	G5126SL		G511SL	G512SL		G5126P	G5126S	
	5/IEC	G5226SL		G521SL	G522SL		G5226P	G5226S	
									

Monitor Manufacturers 监护仪厂家	Leads/label 导联数	Trunk cable 主电缆	Trunk Cable yoke# 分线盒	Leadwire-Pinch 夹式导联线	Leadwire-Snap 扣式导联线	Leadwire Conn.# 导联线接头	Fixed-Pinch 夹式一体化	Fixed-Snap 扣式一体化	ECG Conn. 心电接头
<b>Biolight 宝莱特</b>									
M7000	3/AHA						G3135P	G3135S	CN110
	3/IEC						G3235P	G3235S	
	5/AHA						G5135P	G5135S	
	5/IEC						G5235P	G5235S	
<b>Creative 科瑞康</b>									
PC9000D	5/AHA						G5127P	G5127S	CP6
<b>Comen 科曼</b>									
	3/AHA						G3132P	G3132S	CN101B
	3/IEC						G3232P	G3232S	
<b>General 杰纳瑞</b>									
See Philips	5/AHA						G5124P	G5124S	CN117
<b>M&amp;B 麦邦 Straight</b>									
CD2000	5/AHA						G5128S	G5128S	CN121
									