






Limb Clamps 四肢夹				
P/N 型号	Description 描述	Color 颜色	Package 包装	Photo 图片
K9201	Adult limb clamps electrode with screw, Ag/AgCl, universal ,3mm/4mm, IEC GE OEM: 2029894-001 成人多功能四肢夹, 镀氯化银, 欧标	Red/black/green/yellow 红/黑/绿/黄	4 pcs/pack 4 个/盒	
K9203	Adult limb clamps electrode, Nickel, 3mm/4mm, IEC 成人双功能四肢夹, 镀镍, 欧标	Red/black/green/yellow 红/黑/绿/黄	4 pcs/pack 4 个/盒	
K7201	Pediatric limb clamps electrode with screw, Ag/AgCl, universal ,3mm/4mm, IEC 儿童多功能四肢夹, 镀氯化银, 欧标	Red/black/green/yellow 红/黑/绿/黄	4 pcs/pack 4 个/盒	
Suction Electrode 胸吸球				
P/N 型号	Description 描述	Color 颜色	Package 包装	Photo 图片
K9202	Adult suction electrode with screw, Ag/AgCl, universal, 3mm/4mm Bulb diameter: 26mm GE OEM: 2029892-001 成人多功能胸吸球, 镀氯化银	Blue 蓝色	6 pcs/pack 6 个/盒	
K9204	Adult suction electrode, Nickel, 3mm/4mm Bulb diameter: 26mm GE OEM: 21714401 成人双功能胸吸球, 镀镍	Blue 蓝色	6 pcs/pack 6 个/盒	
K7202	Pediatric suction electrode with screw, Ag/AgCl, universal, 3mm/4mm Bulb diameter: 21mm 儿童多功能胸吸球, 镀氯化银	Blue 蓝色	6 pcs/pack 6 个/盒	