

EKG 心电图机导联线(K)


EKG CABLES&ACCESSORIES 心电图机导联线及配件

- Fixed One-piece, Trunk cable, Leadwire types available
一体化，主电缆，导联线可用
- Latex-free
不含乳胶
- Banana, Needle, Snap, Pinch available
香蕉头，针式，扣式，夹式可用
- AHA and IEC color coding available
美标，欧标可用
- Different types are Compatible with most of patient monitors
多款样式的探头兼容大多数病人监护仪
- Easy to clean and repair
容易清洗和维修



Product types 类型	Photo 图片	P/N Example 型号示例	P/N Guide 型号识别
Fixed one-piece EKG cables 一体化心电图机导联线		<u>K 1 1 13 B</u>	K: EKG cables 1: 10 leads EKG cable 1: AHA (2: IEC) 13: Philips (01: Burdick) B: Banana (N: Needle)
Trunk cable 心电图机主电缆		<u>K 1 1 06 MQ</u>	K: EKG cables 1: 10 leads EKG cable 1: AHA (2: IEC) 06: GE-Marquette MQ: MQ type
Leadwire 心电图机导联线		<u>K 1 1 3 MQ</u>	K: EKG cables 1: 10 leads EKG cable 1: AHA (2: IEC) 3: Banana MQ: MQ type
<p>Warranty: six-month warranty for EKG cables</p> <p>Package: 1pcs/bag</p> <p>保修: 6个月</p> <p>包装: 1个/包</p>			

Monitor Manufacturers 监护仪厂家	Leads/lable 导联数	Fixed Banana 一体化香蕉头	Fixed Needle 一体化插针式	Fixed Snap 一体化扣式	Fixed Pinch 一体化钳式	Trunk cable 主电缆	Leadwire Banana 导联线	Conn. P/N 接头型号
Burdick 伯迪克								
EK-10, Elite, Elite II, E-350, E350i, E550, E560, (1st Generation Eclipse LE/LEi, with serial number starting with 1051), Burdick #007705, Stress cable with detachable leads 007706	10/AHA	K1101B 012-0700-00	K1101N	K1101S	K1101P			CN201
	10/IEC	K1201B	K1201N	K1201S	K1201P			
ATRIA PLUS, ATRIA3000, Eclipse 4 ,Eclipse 400 ,Eclipse 8,Eclipse 800,Eclipse 850 ,Eclipse 850i,Eclipse Plus,Eclipse LE II BURDICK with serial number starting with 1067 or 92304	10/AHA	K1102B 007704 012-0844-01	K1102N	K1102S	K1102P			CN202
	10/IEC	K1202B	K1202N	K1202S	K1202P			
Fukuda Denshi 福田								
FX101 FX102F	10/AHA	K1103B	K1103N	K1103S	K1103P			CN203
	10/IEC	K1203B	K1203N	K1203S	K1203P			
Fukuda Me								
KP-500	10/AHA	K1104B	K1104N	K1104S	K1104P			CN204
	10/IEC	K1204B	K1204N	K1204S	K1204P			
GE-Marquette 马奎								
MAC 500,MAC 1100,MAC 1200,MAC 1200ST	10/AHA	K1106B 2029893-001	K1106N	K1106S	K1106P	K1106MQ 22341809	K113MQ 38401817	CN206
	10/IEC	K1206B 2029890-001	K1206N	K1206S	K1206P	K1206MQ 22341808	K123MQ 38401816	
Kenz 铃谦								
16 male pins, Cardioline (OEM P/N: PC-104)	10/AHA	K1107B	K1107N	K1107S	K1107P			CN207
	10/IEC	K1207B	K1207N	K1207S	K1207P			
(OEM P/N: PC-109)	10/AHA	K1108B	K1108N	K1108S	K1108P			CN208
	10/IEC	K1208B	K1208N	K1208S	K1208P			
Mortara 摩托拉								
	10/AHA	K1107B	K1107N	K1107S	K1107P			CN207
	10/IEC	K1207B	K1207N	K1207S	K1207P			
Nihon Kohden 日本光电								
(Biocare,Dongjiang,Edan, Carewell, Comen, Osen, Spring, Zoncare, Contec) 邦健/东江/理邦/凯沃尔/科曼/奥升/施博瑞/中旗康泰								
ECG-1250,ECG-1350, ECG-9101, ECG-9130, ECG-9132, ECG-9620,PEA-1110 (No resistor)	10/AHA	K1110B	K1110N	K1110S	K1110P			CN210
	10/IEC	K1210B BJ-902D	K1210N BA-901D	K1210S	K1210P BJ-903D			
4.7k resistor	10/AHA	K1111B	K1111N	K1111S	K1111P			CN210
	10/IEC	K1211B	K1211N	K1211S	K1211P			

Monitor Manufacturers 监护仪厂家	Leads/lable 导联数	Fixed Banana 一体化香蕉头	Fixed Needle 一体化插针式	Fixed Snap 一体化扣式	Fixed Pinch 一体化钳式	Trunk cable 主电缆	Leadwire Banana 导联线	Conn. P/N 接头型号
Shanghai Nihon Kohden 上海光电								
	10/AHA	K1112B	K1112N	K1112S	K1112P			CN211
	10/IEC	K1212B	K1212N	K1212S	K1212P			
		BD-615E						
Philips/HP 飞利浦								
M1770A,M1771A,M1772A,M2662A Pagewrite100/200/200i	10/AHA	K1113B	K1113N	K1113S	K1113P			CN212
		M3702C						
	10/IEC	K1213B	K1213N	K1213S	K1213P			
		M3703C						
PageWriter Trim III	10/AHA						K111PH	
	10/IEC						K121PH	
							989803129191	
Schiller 席勒								
AT/1,AT/2,AT/3,AT/4,AT/5,AT/6,AT/10,AT/60,CS6,CS12,C S100-CS200,CV6,P80	10/AHA	K1114B	K1114N	K1114S	K1114P			CN213
	10/IEC	K1214B	K1214N	K1214S	K1214P			
Siemens/Hellige								
	10/AHA	K1115B	K1115N	K1115S	K1115P			CN215
	10/IEC	K1215B	K1215N	K1215S	K1215P			
								

Multi-cable 多功能电缆(M)

Drager-Siemens Compatible MultiMed 5 Pod cable

德尔格-西门子多功能电缆

P/N: M-SM-SGT5

Description: Drager-Siemens MultiMed 5 Pod cable, 2.5m, Compatible with SC6000, 7000, 8000, 9000, DRAGER Infinity, Delta and Delta XL

Parameter: Spo2, ECG(AHA&IEC), Temp

OEM P/N: 3368391 Or 5950196

