



NEBULIZER AND COMPRESSOR 雾化器和压缩雾化吸入机

FROM MEDKE IN SHENZHEN 源自 深圳• 曼克

All ME Nebulizers are designed to be efficient and easy to use with high quality and great price based on the principle of air flow. Its Nebulizers are compatible with most of compressors and air source. Patients are more comfortable because of light weight of material while treatment. Avoid cross-contamination by single patient use.

The TOR (Total Output Rate) is more 350mg/min and particles are small as 3.0um MMD that is more 75%.

ME 型雾化器根据气体流动原理设计，产品效率高，使用简单，性价比高。能够兼容市面上大部分的压缩雾化吸入机及供氧气源。产品轻巧，让病人治疗时更轻松。单病人使用，防止交叉感染。

总输出量(TOR)达到 350 毫克/分钟，总体中位直径(MMD)接近于 3.0 微米，直径 < 5um 的雾粒浓度超过 75%。

ME9531B

Adult Nebulizer kit with mask

[Nebulizer cup, L connector,
Mask, Air tube]

1pcs/bag, 50pcs/carton, 103g/bag

成人面罩雾化器套件

[雾化杯, 弯管接头, 成人面罩, 空气导管]
1 个/包, 50 个/箱, 103 克/包



ME7531G

Pediatric Nebulizer kit with mask

[Nebulizer cup, L connector,
Mask, Air tube]

1pcs/bag, 50pcs/carton, 92g/bag

儿童面罩雾化器套件

[雾化杯, 弯管接头, 儿童面罩, 空气导管]
1 个/包, 50 个/箱, 92 克/包



ME9501B

Adult/Pediatric (>3 years old) Nebulizer kit with mouthpiece





[Nebulizer cup, L connector,
Mouthpiece, Air tube]

1pcs/bag, 50pcs/carton, 75g/bag

成人/儿童(三岁以上)口含雾化器套件

[雾化杯, 弯管接头, 成人/儿童口含管,
空气导管] 1 个/包, 50 个/箱, 75 克/包



Description 描述	Photo 图片	Description 描述	Photo 图片
<p>NP01B Blue Nebulizer cup (medical PP), connect L connector, For mask and mouthpiece</p> <p>蓝色雾化杯(医用 PP 料), 连接弯管接头, 用于面罩和口含管</p>		<p>NP01G Orange Nebulizer cup (medical PP), connect L connector, For mask and mouthpiece</p> <p>橙色雾化杯(医用 PP 料), 连接弯管接头, 用于面罩和口含管</p>	
<p>NP02B Adult mask (medical PVC) with blue Mesh Fabric</p> <p>成人面罩(医用 PVC 料), 带蓝色织带</p>		<p>NP02G Pediatric mask (medical PVC) with orange Mesh Fabric</p> <p>儿童面罩(医用 PVC 料), 带橙色织带</p>	
<p>NP03W Adult/pediatric mouthpiece (medical PP)</p> <p>成人/儿童口含管(医用 PP 料)</p>		<p>NP05W L connector (medical PP), For Adult, Pediatric mask and mouthpiece</p> <p>弯管接头(医用 PP 料), 用于连接成人、儿童面罩及口含管</p>	
<p>NP06B Air tube with blue connectors (medical PVC), 1.2m Connector diameter of air source side: 7.6mm</p> <p>蓝色空气导管(医用 PVC 料), 1.2 米, 带双接头 气源端接口直径: 7.6mm</p>		<p>NP06G Air tube with orange connectors (medical PVC), 1.2m Connector diameter of air source side: 7.6mm</p> <p>橙色空气导管(医用 PVC 料), 1.2 米, 带双接头 气源端接口直径: 7.6mm</p>	



Compressor Nebulizer System

压缩雾化吸入机



- ◆ complete nebulizer kit included
- ◆ Effective Medication Delivery
- ◆ Two storage compartments
- ◆ One-button operation
- ◆ Easy carry handle
- ◆ Ideal for all ages
- ◆ Easy to clean
- ◆ Light weight
- ◆ Latex free

Model: **ME900**
 Net Weight: 1.5kg
 Gross Weight: 1.7kg
 Volume: 33.5cm*7.5cm*21.5cm

Shenzhen MedKe Technology Co., Ltd

A: 4/F, Bldg.A1, Anle Ind. Zone, Hangcheng RD., Baoan, Shenzhen, China 518126
 T: +86 755 23463462 23463463
 F: +86 755 29553084
 E: info@medke.com
 W: www.medke.com

深圳市曼克科技有限公司

地址: 深圳市宝安区固戍航城大道安乐工业区 A1 栋 4 楼 518126
 电话: 0755 23463462 23463463
 传真: 0755-29553084
 邮箱: info@medke.com
 网址: www.medke.com
 淘宝: medke.taobao.com

Ref: ME201404

Distributed by:

