

NIBP 无创血压(C,H)

NIBP CUFFS, CONNECTORS AND ADAPTER HOSES

无创血压，接头及血压连接管

- Reusable and disposable cuffs are available with Single and double tube
重复性，一次性，单双管袖带可用
- Easy to clean
容易清洗
- Latex free
不含乳胶






Product types 产品类型	Photo 图片	P/N Example 型号示例	P/N Guide 型号识别
NIBP disposable cuffs 一次性血压袖带		<u>C 0 1 06</u>	C: NIBP cuffs 0: Disposable cuff 1: single tube (2: two tubes) 06: cuff size
NIBP reusable cuffs with bag 重复性有囊血压袖带		<u>C 1 7 1 1</u>	C: NIBP cuffs 1: Reusable cuff nylon ,6 Reusable cuff PU nylon 7: Adult (6: Pediatric; 5: Infant; 3: Neonate; 8: Large Adult; 9: Thigh) 1: single tube (2: two tubes) 1: with bag
NIBP reusable cuffs without bag 重复性无囊血压袖带		<u>C 1 7 1 0</u>	C: NIBP cuffs 1: Reusable cuff 7: Adult (6: Pediatric; 5: Infant; 3: Neonate; 8: Large Adult; 9: Thigh) 1: single tube (2: two tubes) 0: without bag
NIBP adapter hoses 血压连接管		<u>H 39 07 D</u>	H: NIBP adapter hoses 39 CB39 connector 07: CB07 connector D: Double tubes
Warranty: three-month warranty for NIBP cuffs; six-month warranty for NIBP Adapter Hoses;			
Package: 1pcs/bag			
保修: 重复性袖带 3 个月; 血压连接管 6 个月 包装: 1 个/包			

Disposable NIBP Cuffs 一次性血压袖带

Size 尺寸	Single tube P/N 单管型号	Limb Circumference 臂周长	Double tube P/N 双管型号	Limb Circumference 臂周长	Package 包装
Neonate	C0101-44	3.3-5.6cm	C0201-4444	3.3-5.6cm	1pcs/Bag
Neonate	C0102-44	4.2-7.1cm	C0202-4444	4.2-7.1cm	1pcs/Bag
Neonate	C0103-44	5-10.5cm	C0203-4444	5-10.5cm	1pcs/Bag
Neonate	C0104-44	6.9-11.7cm	C0204-4444	6.9-11.7cm	1pcs/Bag
Neonate	C0105-44	8.9-15cm	C0205-4444	8.9-15cm	1pcs/Bag
Neonate	C0106-44	7.7-10.5cm	C0206-4444	7.7-10.5cm	1pcs/Bag
Infant	C0107-43/C0107-34/C0107-32	9.8-13.3cm	C0207-3636/C0207-3232/C0207-3132	9.8-13.3cm	1pcs/Bag
Small Child	C0108-43/C0108-34/C0108-32	12.4-16.8cm	C0208-3636/C0208-3232/C0208-3132	12.4-16.8cm	1pcs/Bag
Child	C0109-43/C0109-34/C0109-32	15.8-21.3cm	C0209-3636/C0208-3232/C0209-3132	15.8-21.3cm	1pcs/Bag
Small Adult	C0110-43/C0110-34/C0110-32	20-27cm	C0210-3636/C0210-3232/C0210-3132	20-27cm	1pcs/Bag
Adult	C0111-43/C0111-34/C0111-32	25.3-34.3cm	C0211-3636/C0211-3232/C0211-3132	25.3-34.3cm	1pcs/Bag
Large Adult	C0112-43/C0112-34/C0111-32	32.1-43.4cm	C0212-3636/C0212-3232/C0212-3132	32.1-43.4cm	1pcs/Bag
Thigh	C0113-43/C0113-34/C0111-32	46-66cm	C0213-3636/C0213-3232/C0213-3132	46-66cm	1pcs/Bag

Please match right plastic NIBP Connectors for cuffs

Reusable NIBP Cuffs 重复血压袖带

Size 尺寸	Limb Circumference 臂周长	Single tube Blade 单管有囊	Double tube Blade 双管有囊	Photo 图片	Remarks 备注
Adult	25-35cm	C1711	C1721		Please match right NIBP Connectors for cuffs 1pcs/Bag
Pediatric	18-26cm	C1611	C1621		
Infant	10-19cm	C1511	C1521		
Neonate	6-11cm	C1311	C1321		
Large Adult	33-47cm	C1811	C1821		
Thigh	46-66cm	C1911	C1921		
Size 尺寸	Limb Circumference 臂周长	Single tube Bladeless 单管无囊	Double tube Bladeless 双管无囊	Photo 图片	Remarks 备注
Adult	25-35cm	C6711	C6721		Please match right NIBP Connectors for cuffs 1pcs/Bag
Pediatric	18-26cm	C6611	C6621		
Infant	10-19cm	C6511	C6521		
Neonate	6-11cm	C6311	C6321		
Large Adult	33-47cm	C6811	C6821		
Thigh	46-66cm	C6911	C6921		
Size 尺寸	Limb Circumference 臂周长	Single tube Bladeless 单管无囊	Double tube Bladeless 双管无囊	Photo 图片	Remarks 备注
Adult	25.3-34.3cm	C1710	C1720		Please match right NIBP Connectors for cuffs 1pcs/Bag
Pediatric	20-27cm	C1610	C1620		
Infant	12.4-16.8cm	C1510	C1520		
Child	15.8-21.3cm	C1410	C1420		
Neonate	6.9-11.7cm	C1310	C1320		
Large Adult	32.1-43.4cm	C1810	C1820		
Thigh	46-66cm	C1910	C1920		

P/N 型号	Description (Blade) 描述 (有囊)	P/N 型号	Description (Bladeless) 描述 (无囊)	With connector 所带接头
C1711-43	Adult, Nylon, single tube, 25-35cm	C1710-43	Adult, TPU&Nylon, single tube, 25.3-34.3cm	
C1611-43	Pediatric, Nylon, single tube, 18-26cm	C1610-43	Pediatric, TPU&Nylon, single tube, 20-27cm	
C1511-43	Infant, Nylon, single tube, 10-19cm	C1510-43	Infant, TPU&Nylon, single tube, 12.4-16.8cm	
		C1410-43	Child, TPU&Nylon, single tube, 15.8-21.3cm	
C1311-43	Neonate, Nylon, single tube, 6-11cm	C1310-43	Neonate, TPU&Nylon, single tube, 6.9-11.7cm	
C1811-43	Large Adult, Nylon, single tube, 33-47cm	C1810-43	Large Adult, TPU&Nylon, single tube, 32.1-43.4cm	
C1911-43	Thigh, Nylon, single tube, 46-66cm	C1910-43	Thigh, TPU&Nylon, single tube, 46-66cm	
C1711-33	Adult, Nylon, single tube, 25-35cm	C1710-33	Adult, TPU&Nylon, single tube, 25.3-34.3cm	
C1611-33	Pediatric, Nylon, single tube, 18-26cm	C1610-33	Pediatric, TPU&Nylon, single tube, 20-27cm	
C1511-33	Infant, Nylon, single tube, 10-19cm	C1510-33	Infant, TPU&Nylon, single tube, 12.4-16.8cm	
		C1410-33	Child, TPU&Nylon, single tube, 15.8-21.3cm	
C1311-33	Neonate, Nylon, single tube, 6-11cm	C1310-33	Neonate, TPU&Nylon, single tube, 6.9-11.7cm	
C1811-33	Large Adult, Nylon, single tube, 33-47cm	C1810-33	Large Adult, TPU&Nylon, single tube, 32.1-43.4cm	
C1911-33	Thigh, Nylon, single tube, 46-66cm	C1910-33	Thigh, TPU&Nylon, single tube, 46-66cm	
C1711-32	Adult, Nylon, single tube, 25-35cm	C1710-32	Adult, TPU&Nylon, single tube, 25.3-34.3cm	
C1611-32	Pediatric, Nylon, single tube, 18-26cm	C1610-32	Pediatric, TPU&Nylon, single tube, 20-27cm	
C1511-32	Infant, Nylon, single tube, 10-19cm	C1510-32	Infant, TPU&Nylon, single tube, 12.4-16.8cm	
		C1410-32	Child, TPU&Nylon, single tube, 15.8-21.3cm	
C1311-32	Neonate, Nylon, single tube, 6-11cm	C1310-32	Neonate, TPU&Nylon, single tube, 6.9-11.7cm	
C1811-32	Large Adult, Nylon, single tube, 33-47cm	C1810-32	Large Adult, TPU&Nylon, single tube, 32.1-43.4cm	
C1911-32	Thigh, Nylon, single tube, 46-66cm	C1910-32	Thigh, TPU&Nylon, single tube, 46-66cm	
C1721-3636	Adult, Nylon, double tube, 25-35cm	C1720-3636	Adult, TPU&Nylon, double tube, 25.3-34.3cm	
C1621-3636	Pediatric, Nylon, double tube, 18-26cm	C1620-3636	Pediatric, TPU&Nylon, double tube, 20-27cm	
C1521-3636	Infant, Nylon, double tube, 10-19cm	C1520-3636	Infant, TPU&Nylon, double tube, 12.4-16.8cm	
		C1420-3636	Child, TPU&Nylon, double tube, 15.8-21.3cm	
C1321-3636	Neonate, Nylon, double tube, 6-11cm	C1320-3636	Neonate, TPU&Nylon, double tube, 6.9-11.7cm	
C1821-3636	Large Adult, Nylon, double tube, 33-47cm	C1820-3636	Large Adult, TPU&Nylon, double tube, 32.1-43.4cm	
C1921-3636	Thigh, Nylon, double tube, 46-66cm	C1920-3636	Thigh, TPU&Nylon, double tube, 46-66cm	
C1721-3232	Adult, Nylon, double tube, 25-35cm	C1720-3232	Adult, TPU&Nylon, double tube, 25.3-34.3cm	
C1621-3232	Pediatric, Nylon, double tube, 18-26cm	C1620-3232	Pediatric, TPU&Nylon, double tube, 20-27cm	
C1521-3232	Infant, Nylon, double tube, 10-19cm	C1520-3232	Infant, TPU&Nylon, double tube, 12.4-16.8cm	
		C1420-3232	Child, TPU&Nylon, double tube, 15.8-21.3cm	
C1321-3232	Neonate, Nylon, double tube, 6-11cm	C1320-3232	Neonate, TPU&Nylon, double tube, 6.9-11.7cm	
C1821-3232	Large Adult, Nylon, double tube, 33-47cm	C1820-3232	Large Adult, TPU&Nylon, double tube, 32.1-43.4cm	
C1921-3232	Thigh, Nylon, double tube, 46-66cm	C1920-3232	Thigh, TPU&Nylon, double tube, 46-66cm	
C1721-3132	Adult, Nylon, double tube, 25-35cm	C1720-3132	Adult, TPU&Nylon, double tube, 25.3-34.3cm	
C1621-3132	Pediatric, Nylon, double tube, 18-26cm	C1620-3132	Pediatric, TPU&Nylon, double tube, 20-27cm	
C1521-3132	Infant, Nylon, double tube, 10-19cm	C1520-3132	Infant, TPU&Nylon, double tube, 12.4-16.8cm	
		C1420-3132	Child, TPU&Nylon, double tube, 15.8-21.3cm	
C1321-3132	Neonate, Nylon, double tube, 6-11cm	C1320-3132	Neonate, TPU&Nylon, double tube, 6.9-11.7cm	
C1821-3132	Large Adult, Nylon, double tube, 33-47cm	C1820-3132	Large Adult, TPU&Nylon, double tube, 32.1-43.4cm	
C1921-3132	Thigh, Nylon, double tube, 46-66cm	C1920-3132	Thigh, TPU&Nylon, double tube, 46-66cm	