












































































P/N 型号	CB01	CB02	CB03	CB04	CB05	CB06	CB07
Photo 图片							
Match 匹配	CB01->CB02	CB02->CB01	CB03-> CB04	CB04-> CB03	CB05-> CB06	CB06-> CB05	CB07-> CB36/CB16
Use with 适用于	For Neonate air hose conn.	CAS, Drager, Kontron, Physio control, Schiller /Bruker/Odam	Mennen (Old), Schiller/Bruker/O dam, Spacelabs, Nihon Kohden Monitor plug	BCI, Datascope(Old)	Monitor plug	Colin, Gambro, Tyco(NPB 5500)	Critikon, Invivo, Kontron, Nihon Kohden GE: 300619
P/N 型号	CB08	CB09	CB10	CB11	CB12	CB13	CB14
Photo 图片							
Match 匹配			CB10-> CB11	CB11-> CB10			
Use with 适用于	ID Barbed to ID Barbed, Metal	Philips/HP Monitor plug	Mindray, Schiller GE: 330060	Philips, Siemens, Spacelabs, CSI, Welch Allyn, Datascope GE: 330059 WA: 5082-179	Drager Monitor plug	Monitor plug	PACE TECH Monitor plug
P/N 型号	CB16	CB30	CB31	CB32	CB33	CB34	CB36
Photo 图片							
Match 匹配	CB16-> CB07		CB31-> CB32	CB32->CB31	CB33-> CB34	CB34-> CB33	CB36-> CB07
Use with 适用于	Critikon, Invivo, Kontron, Mennen, Nihon Kohden, Welch Allyn GE: 300664 WA: 5082-164	ID Barbed to ID Barbed, Plastic GE: 300670	Welch Allyn, Marquette, Datex-Ohmeda, Nihon Kohden GE: 330091 WA: 5082-178	Welch Allyn Marquette, Datex, Zoll GE: 330064 WA: 5082-172	BCI, Datascope GE: 300665 WA: 5082-165	BCI, Datascope GE: 300668 WA: 5082-168	Critikon, Invivo, Kontron, Mennen, Nihon Kohden, Welch Allyn GE: 300664 WA: 5082-164
P/N 型号	CB37	CB38	CB39	CB40	CB41	CB42	CB43
Photo 图片							
Match 匹配							CB43-> CB10
Use with 适用于	GE AS/3, CS/3, S/5, Eagle4000, Cardiacap 5 Monitor plug	GE-Dinamap Monitor plug	GE-Dash, Eagle, Solar Monitor plug	Nihon Kohden For adult/pediatric Monitor plug	Nihon Kohden For neonate/infant Monitor plug	for neonate double tube disposable BP cuffs	Philips, Siemens, Spacelabs, CSI, Datascope WA: 5082-184
P/N 型号	CB44	CB45	CB46				
Photo 图片							
Match 匹配							
Use with 适用于	Neonate disposable cuff WA: 5082-176	Use with all NIBP connector					

NIBP ADAPTER HOSES 血压连接管

Monitor Manufacturers/ Model# 监护仪厂家/监护仪型号	Tube 单双管	Monitor side Conn. 监护仪端接头	OEM P/N Adult&Pedia. 连接管原厂型号	P/N for Adult&Pedia. 成人/儿童型号	Cuff side Conn. 袖带端接头	OEM P/N Infant&Neo. 连接管原厂型号	P/N for Infant&Neo. 婴儿/新生儿型号	Cuff side Conn. 袖带端接头
Air shields								
Athena	2	 CB11		H1107D	 CB07		H1142ND	 CB42
Colin 科林								
8800	1	 CB05		H0506S	 CB06		H0501NS	 CB01
CSI								
507,1100	2	 CB37	558	H3707D	 CB07		H3742ND	 CB42
	1	 CB05	709	H0510S	 CB10		H0501NS	 CB01
	1	 CB10	706	H1010S	 CB10		H1001NS	 CB01
507E,508,Scholar II (See Colin)	1	 CB05		H0506S	 CB06		H0501NS	 CB01
Drager-Siemens 德尔格-西门子								
SC6000-7000/8000/9000Infinity	1	 CB12	1275275 EH50U	H1210S	 CB10	2870298	H1201NS	 CB01
GE-Critikon 强生								
Dinamap 1846SX,Dinamap 8100,Dinamap 9300XL	2	 CB07	008841 E9050JY	H0707D	 CB07	008840	H0742ND	 CB42
(See CSI)	2	 CB37	9461-214	H3707D	 CB07		H3742ND	 CB42

























NIBP ADAPTER HOSES 血压连接管

Monitor Manufacturers/ Model# 监护仪厂家/监护仪型号	Tube 单双管	Monitor side Conn. 监护仪端接头	OEM P/N Adult&Pedia. 连接管原厂型号	P/N for Adult&Pedia. 成人/儿童型号	Cuff side Conn. 袖带端接头	OEM P/N Infant&Neo. 连接管原厂型号	P/N for Infant&Neo. 婴儿/新生儿型号	Cuff side Conn. 袖带端接头
GE-Datex Ohmeda 德恩								
AS/3,CS/3,Cardiacap 5,Cardiacap II, Light Monitor S/5	2		877235	H37312D	 	877514	H3742ND	
GE-Marquette 马奎								
Eagle 4000, Tram 100/200/250/300/350, Corometrics 120	2		9461-203 E9001TE	H3907D		877514	H3942ND	
Eagle 3000, DASH 2000/3000/4000, Tram 400/450/451, Tram 800/850/851, Corometrics 120	2		414873-00 1 E9001TC 2017008-0 01	H3931D		414874-0 01 2017009- 001	H3942ND	
Dinamap MPS Series, Dinamap PRO 100/110/200/210/300/310/400 Dinamap ProCare 100/120/200/220/300/320/400	2		107365	H3807D		107368	H3842ND	
Invivo / MDE 英威								
1400,5600 (See Nihon kohden)	2			H0707D			H0742ND	
Kontron (AVCO) 康强								
Minimon 7137B, Minimon 7137Plus (See Critikon)	2			H0707D			H0742ND	
Medtronic-Physio Contro 美敦力								
Lifepak 12	1			H1110S			H1101NS	
Mindray 迈瑞								
PM-600, PM-6000, PM-7000, PM-8 000, PM-8000 Express, PM-9000 Express, VS-800 (See CSI)	1		M3918A	H1010S			H1001NS	

NIBP ADAPTER HOSES 血压连接管

Monitor Manufacturers/ Model# 监护仪厂家/监护仪型号	Tube 单双管	Monitor side Conn. 监护仪端接头	OEM P/N Adult&Pedia. 连接管原厂型号	P/N for Adult&Pedia. 成人/儿童型号	Cuff side Conn. 袖带端接头	OEM P/N Infant&Neo. 连接管原厂型号	P/N for Infant&Neo. 婴儿/新生儿型号	Cuff side Conn. 袖带端接头
Mindray 迈瑞								
BeneView T5,BeneView T6,BeneView T8,MPM Module (See Philips)	1			H0910S			H0901NS	
Mindray-Datascope								
3000,Accutorr 1,Accutorr 2, Accutorr 3,Accutorr 4,Accutorr Plus,ExpertPassport,Passport 2, Passport 2LT,Spectrum (See CSI)	1		0683-04-006	H0510S			H0501NS	
2100,2200I,Accutorr 1,Accutorr 2, Accutorr 3,Accutorr 4,Accutorr Plus,Passort 3L,Passport,Passport 2,Passport 2LT,Passport 5L,Spectrum	1		0683-04-007	H0610S			H0601NS	
AccuTorr Plus,Duo,Passport 2,Passport 2LT,Spectrum,Spectrum OR Trio	1		0683-04-003 0683-04-004	H1010S			H1001NS	
Nellcor								
3100,4000	1			H1010S			H1001NS	
Nihon Kohden 光电								
AP-860PA Module, BSM-2101,BSM-2102,BSM-8301, BSM-8302,BSM-8500,BSM-8502, BSM-8800,Lifescop 14, Lifescop 9 Lifescop L	2		6143-001321	H0707D			H0742ND	
AP-851PA Module,BSM-8300, BSM-8301,BSM-8302,Lifescop 9	1		110088	H0304S			H0301NS	
CB40 is different from CB41	1		YN-900P YN-901P	H4032D		YN-920P YN-921P	H4101D	
Philips/HP 飞利浦								
78352A,78352C,78354A,78354C M1008A,M1008B,M3000	1		M1599B M1598B	H0910S		M1597B M1596B	H0901NS	

NIBP ADAPTER HOSES 血压连接管

Monitor Manufacturers/ Model# 监护仪厂家/监护仪型号	Tube 单双管	Monitor side Conn. 监护仪接头	OEM P/N Adult&Pedia. 连接管原厂型号	P/N for Adult&Pedia. 成人/儿童型号	Cuff side Conn. 袖带接头	OEM P/N Infant&Neo. 连接管原厂型号	P/N for Infant&Neo. 婴儿/新生儿型号	Cuff side Conn. 袖带接头
Philips/HP 飞利浦								
Compact A1, Compact A3	1	 CB10	M3918A	H1010S	 CB10		H1001NS	 CB01
Spacelabs Medical 太空								
Ultraview 90496, Ultraview 91496	1	 CB10	714-0018-00	H1010S	 CB10		H1001NS	 CB01
Ultraview 90496, Ultraview 91496	1	 CB11	714-0019-01	H1134S	 CB34		H1101NS	 CB01
All (dual tube) (See Nihon Kohden)	2	 CB07		H0707D	 CB07		H0742ND	 CB42
All (single luer) except PCI (See Datascope)	1	 CB06		H0610S	 CB10		H0601NS	 CB01
Welch Allyn 伟伦								
Atlas, Protocol (See Nihon Kohden)	1	 CB03		H0304S	 CB04		H0301NS	 CB01
Encore 202EL, Encore 204EL, Encore 206EL, Propaq 102EL, Propaq 104EL, Propaq 106EL, Protocol	1	 CB07		H0731S	 CB31		H0701NS	 CB01
Zoll Medical								
Masim Series (See Datascope)	1	 CB06		H0610S	 CB10		H0601NS	 CB01