











- Finger clip, Soft tip, Wrap, Multi-Y, Ear clip types available  
指夹，指套，包裹，多功能 Y 型，耳夹可用
- Latex-free  
不含乳胶
- Available in adult, pediatric, infant and neonate sizes  
成人，儿童，婴儿，新生儿可用
- Different types are Compatible with most of patient monitors  
多款样式的探头兼容大多数病人监护仪
- High performance/price ratio for Reusable SpO2 sensor to save cost  
高性价比的重复性血氧探头，节约成本
- Easy to clean and repair  
容易清洗和维修





Product types 产品类型	Photo 图片	P/N Example 型号示例	P/N Guide 型号识别
Reusable SpO2 Sensors 重复性血氧探头		<u>P 8 3 25A</u>	<b>P:</b> SpO2 <b>8:</b> Adult soft (9: Adult clip; 7: Pediatric clip; 6: Pediatric soft; 5: Neonate wrap; 4: Multi Y; 3: Ear clip) <b>3:</b> 3.0m (1: 1.0m) <b>25A:</b> Philips 8 pins (25: Philips 12 pins; 19: Nellcor 7 pins; 19A: Nellcor 9 pins)
	      		
SpO2 Adapter Cables 血氧转接线		<u>P 0 2 21B</u>	<b>P:</b> SpO2 <b>0:</b> SpO2 Adapter Cables <b>2:</b> 2.2m (1: 1.0m; 3: 3.0m) <b>21B:</b> NIHON KOHDEN
<b>Warranty: one-year warranty for reusable sensors and cables</b> <b>Package: 1pcs/bag</b> 保修：一年 包装：1 个/包			

Monitor Manufacturers 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能 Y 型	Ear clip 成人耳夹	Len. 长度	Conn. 接头	Adapter Cable 转接线
<b>Armstrong Medical</b>										
AD-900, 1000, 2000, 3000, 4000(BCI tech)	P9303A	P8303A	P7303A	P6303A	P5303A	P4303A	P3303A	3m	CN81	Direct 直连
<b>Artema 阿体玛</b>										
Diascope	P9301	P8301	P7301	P6301	P5301	P4301	P3301	3m	CN01	Direct 直连
S&W	P9301A	P8301A	P7301A	P6301A	P5301A	P4301A	P3301A	3m	CN02	Direct 直连
	SL033057									
<b>Baxter/Simed 百特</b>										
ASAT	P9302	P8302	P7302	P6302	P5302	P4302	P3302	3m	CN03	Direct 直连
	600015									
Monarch/Explorer	P9319	P8319	P7319	P6319	P5319	P4319	P3319	3m	CN30	Direct 直连
S-50, S100, S-100e(Simed)	P9302	P8302	P7302	P6302	P5302	P4302	P3302	3m	CN03	Direct 直连
<b>BCI/Smiths</b>										
3100, 3101, 6100, 9100, Advisor, Autocorr, Mini Torr, Mini Torr Plus	P9303	P8303	P7303	P6303	P5303	P4303	P3303	3m	CP7	Direct 直连
3300, 3301, 3302, 3303, 3304, 3401, 3180	P9303A	P8303A	P7303A	P6303A	P5303A	P4303A	P3303A	3m	CN81	Direct 直连
	3044/3444				3025/3026	3043				
	P9103A	P8103A	P7103A	P6103A	P5103A	P4103A	P3103A	1m	CN81	See cable
<b>Biosys</b>										
BPM-200, Vital Sign Monitor SENTRY	P9304	P8304	P7304	P6304	P5304	P4304	P3304	3m	CN05	Direct 直连
<b>Bird</b>										
4400	P9307	P8307	P7307	P6307	P5307	P4307	P3307	3m	CM5	Direct 直连
<b>Bionet</b>										
	P9306	P8306	P7306	P6306	P5306	P4306	P3306	3m	CP7-1	Direct 直连
	P9306A	P8306A	P7306A	P6306A	P5306A	P4306A	P3306A	3m	CN63	Direct 直连
	P9306B	P8306B	P7306B	P6306B	P5306B	P4306B	P3306B	3m	CN78	Direct 直连
<b>CAS Medical</b>										
OscilloMate 9002, 9303, 9302S(Nonin tech)	P9322B	P8322B	P7322B	P6322B	P5322B	P4322B	P3322B	3m	CN83	Direct 直连
740-750(Masim tech)	P9315	P8315	P7315	P6315	P5315	P4315	P3315	3m	CN22	Direct 直连
<b>Charmcare</b>										
CX100(Nellcor tech)	P9319	P8319	P7319	P6319	P5319	P4319	P3319	3m	CN30	Direct 直连
<b>Colin 科林</b>										
BP88/306/508/608/BX10ne(Nellcor tech)	P9319	P8319	P7319	P6319	P5319	P4319	P3319	3m	CN30	Direct 直连
BP 88S	P9319B	P8319B	P7319B	P6319B	P5319B	P4319B	P3319B	3m	CN22	Direct 直连
BX10ma(Masim tech)	P9315	P8315	P7315	P6315	P5315	P4315	P3315	3m	CN22	Direct 直连
<b>Corpuls</b>										
08/16(BCI tech)	P9303A	P8303A	P7303A	P6303A	P5303A	P4303A	P3303A	3m	CN81	Direct 直连
08/16(M)(Nellcor tech)	P9319	P8319	P7319	P6319	P5319	P4319	P3319	3m	CN30	Direct 直连




## REUSABLE SPO2 SENSORS 重复性血氧探头

Monitor Manufacturers 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能 Y 型	Ear clip 成人耳夹	Len. 长度	Conn. 接头	Adapter Cable 转接线
<b>CSI</b>										
Straight connector	P9107B	P8107B	P7107B	P6107B	P5107B	P4107B	P3107B	1m	CN47	Direct 直连
501, 503, 504, 507, 508, 581, 1100, 8100,Poet (602-3A), Scholar	P9307	P8307	P7307	P6307	P5307	P4307	P3307	3m	CM5	Direct 直连
	934-10LN	975AD-10L	1563-10L	975PD-10L		940-10L	1567-10L			
503 Spot, 504Plus, 504DX, 506DXN	P9307A	P8307A	P7307A	P6307A	P5307A	P4307A	P3307A	3m	CN06	Direct 直连
90 degree angle	934-10DN	975AD-10	1563-10D	975PD-10		940-10D	1567-10D			
90 degree angle	P9107A	P8107A	P7107A	P6107A	P5107A	P4107A	P3107A	1m	CN06	See cable
	934SDN	975AD-3	1563-SD	975PD-3		940SD	1567SD			
<b>Dolphin 海豚</b>										
2100, 2150 (Digital)	P9162A	P8162A	P7162A	P6162A	P5162A	P4162A	P3162A	1m	CN71	See cable
	P9362	P8362	P7362	P6362	P5362	P4362	P3362	3m	CN22	Direct 直连
<b>Dräger-Siemens 德尔格-西门子</b>										
Dräger Oxipac 2000, 2500(BCI tech)	P9303A	P8303A	P7303A	P6303A	P5303A	P4303A	P3303A	3m	CN81	Direct 直连
Vitalert 1000, Cicero, GammaXL, DeltaXL, Kappa	P9309	P8309	P7309	P6309	P5309	P4309	P3309	3m	CM8	Direct 直连
Sirecust, Sirem Spo2	P9309A	P8309A	P7309A	P6309A	P5309A	P4309A	P3309A	3m	CN09	Direct 直连
SC6000, 7000,9000	P9309B	P8309B	P7309B	P6309B	P5309B	P4309B	P3309B	3m	CN10	Direct 直连
<b>GE-Corometrics 科洛</b>										
Corometrics 188,129,129F,128, 128F	P9319	P8319	P7319	P6319	P5319	P4319	P3319	3m	CN30	Direct 直连
Corometrics 120, 510, 511, 118,556	P9310	P8310	P7310	P6310	P5310	P4310	P3310	3m	CN11	Direct 直连
<b>GE-Critikon 强生</b>										
Dinamap 8700/8710/8720 Plus, Oxysuttle	P9310A	P8310A	P7310A	P6310A	P5310A	P4310A	P3310A	3m	CM7	Direct 直连
Dinamap 9700 / 9710 / 9720, Plus,Compact, MPS, Select	P9319	P8319	P7319	P6319	P5319	P4319	P3319	3m	CN30	Direct 直连
Dinamap ProSeries V1 and Pro 1000 V2 (Nelcor tech)	P9319B	P8319B	P7319B	P6319B	P5319B	P4319B	P3319B	3m	CN22	Direct 直连
ProCare; ProSeries V2; Pro 1000 V3(Masim)	P9315	P8315	P7315	P6315	P5315	P4315	P3315	3m	CN22	Direct 直连
Same as above(oximax)	P9119A	P8119A	P7119A	P6119A	P5119A	P4119A	P3119A	1m	CN31	P0319C
<b>GE-Datex Ohmeda 德恩</b>										
S/5, AS/3, CS/3 Modular Monitors; with M-series modules; Cardiocap/5, Cardiocap	P9310B	P8310B	P7310B	P6310B	P5310B	P4310B	P3310B	3m	CN12	Direct 直连
	OXY-F4-N						OXY-E4-N			
S/5 Light Monitor;	P9110C	P8110C	P7110C	P6110C	P5110C	P4110C	P3110C	1m	CN04	P0210B
Capnomac Ultima Satllite Plus II	OXY-F-DB				OXY-W-DB		OXY-E-DB			
M-NSAT Modules for AS/3, CS/3,S/5 Modular	P9319B	P8319B	P7319B	P6319B	P5319B	P4319B	P3319B	3m	CN22	Direct 直连
E-NSAT Modules for i4 Modular Monitors	P9119A	P8119A	P7119A	P6119A	P5119A	P4119A	P3119A	1m	CN31	P0319C


Monitor Manufacturers 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能 Y 型	Ear clip 成人耳夹	Len. 长度	Conn. 接头	Adapter Cable 转接线
<b>GE-Marquette 马奎</b>										
Dash 2000/3000/4000(Nellcor); Solar 8000/8000M,i/9500 with Tram Module 451/451N/851/851N; CardioLab /Mac-Lab with Tram Module 451/451N; Corometrics 120/2120 series; Eagle 1000/3000/4000; Responder 3000 (Nellcor tech)	P9310E	P8310E	P7310E	P6310E	P5310E	P4310E	P3310E	3m	CN14	Direct 直连
Dash 3000/4000/5000(OxiMax); Solar 8000/8000M, i9500 with Tram module 451N5/851N5; Dash 3000/4000/5000(OxiMax); Solar 8000/8000M, i9500 with Tram module 451N5/851N5	P9310F	P8310F	P7310F	P6310F	P5310F	P4310F	P3310F	3m	CN15	Direct 直连
Dash 3000/4000/5000 (Masim) Solar 8000/8000M, i/9500 with Tram Module 451M/851M CardioLab /Mac-Lab with Tram Module 451M; Corometrics 259(Masim)	P9310G	P8310G	P7310G	P6310G	P5310G	P4310G	P3310G	3m	CN15	Direct 直连
<b>GE-Hellige 海力格</b>										
Eagle 1000, 3000, SMU, SMK, SMV, cardioserv	P9310D	P8310D	P7310D	P6310D	P5310D	P4310D	P3310D	3m	CN13	Direct 直连
<b>GE-Ohmeda 欧美达</b>										
Ohmeda TuffSat;3775, 37xx Series, 3800, 3900,3900P, 4700 and 5250 respiratory gas monitors;Cardiocap/5 with N-XOSAT option;S/5, AS/3, CS/3 Modular Monitors with M-OSAT module	P9310H	P8310H	P7310H	P6310H	P5310H	P4310H	P3310H	3m	CN16	Direct 直连
	OXY-F4-H				OXY-W4-H		OXY-E4-H			
Ohmeda TuffSat	P9110J	P8110J	P7110J	P6110J	P5110J	P4110J	P3110J	1m	CN04	P0210H
	P9110H	P8110H	P7110H	P6110H	P5110H	P4110H	P3110H	1m	CN16	Direct 直连
	OXY-F1-H									
S/5, i4 Modular Monitors with E-series modules; Corometrics 259 with TruSignal(Ohmeda tech) B20/30/40 Series	P9310K	P8310K	P7310K	P6310K	P5310K	P4310K	P3310K	3m	CN15	Direct 直连
	OXY-F4-GE									
GE-Ohmeda OxyTip sensor	P9110L	P8110L	P7110L	P6110L	P5110L	P4110L	P3110L	1m	CN17	See
	OXY-F-UN				OXY-W-UN	OXY-SE-3				GE cable
Ohmeda TruSat, Ohmeda TruSat TD	P9310M	P8310M	P7310M	P6310M	P5310M	P4310M	P3310M	3m	CN18	Direct 直连
	OXY-F4-MC						OXY-E4-MC			
<b>GE TruSignal</b>										
TruSignal sensor 	P9110S	P8110S	P7110S	P6110S	P5110S	P4110S		1m	CN64	See cable
	TS-F-D	TS-SA-D		TS-SP-D	TS-W-D	TS-SE-3				
							P3110T	1m	CN74	See cable
							TS-E-D			
										



Monitor Manufacturers 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能 Y 型	Ear clip 成人耳夹	Len. 长度	Conn. 接头	Adapter Cable 转接线
<b>Healthdyne</b>										
All	P9322B	P8322B	P7322B	P6322B	P5322B	P4322B	P3322B	3m	CN04	Direct 直连
<b>Infinium</b>										
OMNI 2, OMNI 3; 武汉思创 (Digital)	P9365	P8365	P7365	P6365	P5365	P4365	P3365	3m	CP5R	Direct 直连
OMNI(BCI Tech)	P9365A	P8365A	P7365A	P6365A	P5365A	P4365A	P3365A	3m	CP6	Direct 直连
<b>Innomed</b>										
Oxycard, Innocare-T0	P9322B	P8322B	P7322B	P6322B	P5322B	P4322B	P3322B	3m	CN04	Direct 直连
<b>Invacare</b>										
VitalGaurd	P9303A	P8303A	P7303A	P6303A	P5303A	P4303A	P3303A	3m	CN04	Direct 直连
<b>Invivo 英威</b>										
4500, 4500 Plus, Millenia, Omega 1445, Omni-Trak TVS, Sat Scout, 4500 Plus 2	P9312	P8312	P7312	P6312	P5312	P4312	P3312	3m	CM7	Direct 直连
New generation(Masim)	P9312A	P8312A	P7312A	P6312A	P5312A	P4312A	P3312A	3m	CM7	Direct 直连
<b>Ivy Biomedical</b>										
405P, 405T	P9319	P8319	P7319	P6319	P5319	P4319	P3319	3m	CN30	Direct 直连
Vital-Guard 450C	P9315	P8315	P7315	P6315	P5315	P4315	P3315	3m	CN22	Direct 直连
<b>Kontron 康强</b>										
7840, 7845, Minimon 7137, Micromon, Supermon, Kolormon, Trakmon	P9313 0608010	P8313	P7313	P6313	P5313	P4313	P3313	3m	CN20	Direct 直连
7 pins 80 degree	P9313A	P8313A	P7313A	P6313A	P5313A	P4313A	P3313A	3m	CP7T	Direct 直连
<b>Lohmerier 隆梅尔</b>										
	P9314	P8314	P7314	P6314	P5314	P4314	P3314	3m	CN21	Direct 直连
<b>L&amp;T Medical</b>										
Stellar P, Star 50, Galaxy, Planet, Planet 50, Cosmic	P9303A	P8303A	P7303A	P6303A	P5303A	P4303A	P3303A	3m	CN04	Direct 直连
	P9376A	P8376A	P7376A	P6376A	P5376A	P4376A	P3376A	3m	CN79	Direct 直连
<b>Lynn Medical</b>										
All	P9303A	P8303A	P7303A	P6303A	P5303A	P4303A	P3303A	3m	CN04	Direct 直连
<b>Masim 玛思莫</b>										
	P9315	P8315	P7315	P6315	P5315	P4315	P3315	3m	CN22	Direct 直连
	1396/1969									
	P9115A	P8115A	P7115A	P6115A	P5115A	P4115A	P3115A	1m	CN23	See Masim cable
	1863/LNCS DCI									
	P9115B	P8115B	P7115B	P6115B	P5115B	P4115B	P3115B	1m	CN24	
	LNOP DCI		LNOP DC-IP			1544	LNOP TC-I			
<b>M-LNCS™ DCI</b>	P9115H	P8115H	P7115H	P6115H	P5115H	P4115H	P3115H	1m	CN86	
DB9PIN(M* chip)	P9115F	P8115F	P7115F	P6115F	P5115F	P4115F	P3115F	1m	CN23	
Radical-7,Rad-87,Rad-57,Pronto	P9315G	P8315G	P7315G	P6315G	P5315G	P4315G	P3315G	3m	CN80	Direct 直连


Monitor Manufacturers 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能 Y 型	Ear clip 成人耳夹	Len. 长度	Conn. 接头	Adapter Cable 转接线
<b>M.D.E.</b>										
Escort II OPT30; Prism	P9319	P8319	P7319	P6319	P5319	P4319	P3319	3m	CN30	Direct 直连
<b>Mediana</b>										
YM-2002 Series	P9119A	P8119A	P7119A	P6119A	P5119A	P4119A	P3119A	1m	CN31	P0319C
<b>Medlab</b>										
NANOX2, NANOXeco	P9303A	P8303A	P7303A	P6303A	P5303A	P4303A	P3303A	3m	CN04	Direct 直连
PEARL10, NANOX10, Nanox10C	P9319	P8319	P7319	P6319	P5319	P4319	P3319	3m	CN30	Direct 直连
CAPNOX, VITRO, élan, POX 10, P-OX100L	P9303	P8303	P7303	P6303	P5303	P4303	P3303	3m	CP7	Direct 直连
<b>Medsource</b>										
BPM-200	P9319	P8319	P7319	P6319	P5319	P4319	P3319	3m	CN30	Direct 直连
<b>MEK</b>										
MP 100,110,400,1000	P9316	P8316	P7316	P6316	P5316	P4316	P3316	3m	CN30	Direct 直连
	P9116	P8116	P7116	P6116	P5116	P4116	P3116	1m	CN30	P0219
Normal 8 pins	P9316A	P8316A	P7316A	P6316A	P5316A	P4316A	P3316A	3m	CN65	Direct 直连
	P9316B	P8316B	P7316B	P6316B	P5316B	P4316B	P3316B	3m	CN65	Direct 直连
	P9316C	P8316C	P7316C	P6316C	P5316C	P4316C	P3316C	3m	CN66	Direct 直连
<b>Mennen 美能</b>										
Horizon XL, Mercury 1200,MR 1300, Envoy	P9317	P8317	P7317	P6317	P5317	P4317	P3317	3m	CM10	Direct 直连
Normal 8 pins	P9317A	P8317A	P7317A	P6317A	P5317A	P4317A	P3317A	3m	CN25	Direct 直连
(Masim tech)	P9317B	P8317B	P7317B	P6317B	P5317B	P4317B	P3317B	3m	CN26	Direct 直连
<b>Mindray 迈瑞</b>										
12 Male pins	P9318	P8318	P7318	P6318	P5318	P4318	P3318	3m	CN39	Direct 直连
12 Male pins (Masim tech)	P9318A	P8318A	P7318A	P6318A	P5318A	P4318A	P3318A	3m	CN39	Direct 直连
503,507,508,509 (Criticare tech)	P9318B	P8318B	P7318B	P6318B	P5318B	P4318B	P3318B	3m	CM5	Direct 直连
Normal 5 pins	P9318C	P8318C	P7318C	P6318C	P5318C	P4318C	P3318C	3m	CP5	Direct 直连
Normal 6 pins	P9318D	P8318D	P7318D	P6318D	P5318D	P4318D	P3318D	3m	CP6	Direct 直连
PM7000, 8000, 9000, MEC-1000, MEC-2000(6 pins, 40 degree) 科瑞康 PC-900 PC-9000 M8B	P9318E	P8318E	P7318E	P6318E	P5318E	P4318E	P3318E	3m	CP6R	Direct 直连
6 male pins, 40 degree(Masim tech)	P9318F	P8318F	P7318F	P6318F	P5318F	P4318F	P3318F	3m	CP6R	Direct 直连
6 male pins, 60 degree(Masim tech)	P9318G	P8318G	P7318G	P6318G	P5318G	P4318G	P3318G	3m	CP6S	Direct 直连
T5/T8 (Mindray tech)	P9318H	P8318H	P7318H	P6318H	P5318H	P4318H	P3318H	3m	CN27	Direct 直连
T5/T8(Masim tech)	P9318J	P8318J	P7318J	P6318J	P5318J	P4318J	P3318J	3m	CN28	Direct 直连
T5/T8(Nellcor Oximax)	P9318K	P8318K	P7318K	P6318K	P5318K	P4318K	P3318K	3m	CN29	Direct 直连
	P9318L	P8318L	P7318L	P6318L	P5318L	P4318L	P3318L	3m	CP6T	Direct 直连



Monitor Manufacturers 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能 Y 型	Ear clip 成人耳夹	Len. 长度	Conn. 接头	Adapter Cable 转接线
<b>Mindray-Datascope</b>										
Accusat, Accutor 3/4/SAT/PLUS	P9308	P8308	P7308	P6308	P5308	P4308	P3308	3m	CN07	Direct 直连
Passport(Datascope Tech)	P9108A	P8108A	P7108A	P6108A	P5108A	P4108A	P3108A	1m	CN08	PO208
Passport II/IILT, XG, Accutor Plus, Expert (Nellcor) (Trio/Spectrum w/opt. Nellcor)	P9319B	P8319B	P7319B	P6319B	P5319B	P4319B	P3319B	3m	CN22	Direct 直连
Passport II/IILT, XG, Accutor Plus, Trio,Spectrum (Masim)	P9315	P8315	P7315	P6315	P5315	P4315	P3315	3m	CN22	Direct 直连
	P9115B	P8115B	P7115B	P6115B	P5115B	P4115B	P3115B	1m	CN24	PO215B
<b>Nellcor</b>										
NPB 40, NPB20, NPB75, NPB 180, NPB 185,NPB-190, NPB-195,NPB 3900, NPB 4000	P9319	P8319	P7319	P6319	P5319	P4319	P3319	3m	CN30	Direct 直连
	P9119	P8119	P7119	P6119	P5119	P4119	P3119	1m	CN30	See cable
N550,N560,N595,N600,N751, N5500	P9319A	P8319A	P7319A	P6319A	P5319A	P4319A	P3319A	3m	CN31	Direct 直连
DASH 2500, SOLAR 8000, 8000M, N65 Except CARESCAPE V100	P9119A	P8119A	P7119A	P6119A	P5119A	P4119A	P3119A	1m	CN31	See SpO2 cable
	DS-100A					D-YS				
NPB3000-3100, Symphony, NPB290,295,395	P9319B	P8319B	P7319B	P6319B	P5319B	P4319B	P3319B	3m	CN22	Direct 直连
DOC-10 direct sensor	P9319C	P8319C	P7319C	P6319C	P5319C	P4319C	P3319C	3m	CN46	Direct 直连
<b>Nihon Kohden 光电</b>										
AL-800PA Module; BSM-1101, 2101/2, 3101, 7100, 7200; OLV-1100A, Lifescope	P9321	P8321	P7321	P6321	P5321	P4321	P3321	3m	CN32	Direct 直连
	TL-101S									
BSM-2301K, BSM-2303K, BSM-2351, BSM-2353, BSM-4103, BSM-4113,	P9121A	P8121A	P7121A	P6121A	P5121A	P4121A	P3121A	1m	CN33	PO221B
	TL-201T				TL-220T					
	P9321B	P8321B	P7321B	P6321B	P5321B	P4321B	P3321B	3m	CN34	Direct 直连
Lifescope OPV 1500 1500k	P9321C	P8321C	P7321C	P6321C	P5321C	P4321C	P3321C	3m	CN84	
BSM Series ,Oximax	P9321D	P8321D	P7321D	P6321D	P5321D	P4321D	P3321D	3m	CN85	
<b>Nonin 燕牌</b>										
8604	P9322	P8322	P7322	P6322	P5322	P4322	P3322	3m	CM6	Direct 直连
8500, PamSAT 2500, Avant 9600, Avant 9700, Avant 2120, 9840	P9122B	P8122B	P7122B	P6122B	P5122B	P4122B	P3122B	1m	CN83	See cable
	8000AA-1+8500VI									
<b>Novamatrix 诺医</b>										
505, 510, 511, 515B, 515C, 520, 7000, 7100, Oxypleth	P9323	P8323	P7323	P6323	P5323	P4323	P3323	3m	CN36	Direct 直连
	8776-00									
512/513/710/715,MARS	P9123A	P8123A	P7123A	P6123A	P5123A	P4123A	P3123A	1m	CN37	Direct 直连
	P9123B	P8123B	P7123B	P6123B	P5123B	P4123B	P3123B	1m	CN04	PO223
<b>Pace Tech/Generra</b>										
Vitalmax 800 Series	P9311C	P8311C	P7311C	P6311C	P5311C	P4311C	P3311C	3m	CN19	Direct 直连
<b>Palco 派克</b>										
All	P9324	P8324	P7324	P6324	P5324	P4324	P3324	3m	CN38	Direct 直连
	P9324A	P8324A	P7324A	P6324A	P5324A	P4324A	P3324A	3m	CN67	Direct 直连



Monitor Manufacturers 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能 Y 型	Ear clip 成人耳夹	Len. 长度	Conn. 接头	Adapter Cable 转接线
<b>Philips-HP- Agilent 飞利浦</b>										
78352C, 78354C, 78834C, M1020A, M1025B,M1025A, M1350B,M1722A/B, M1732A/B, M2475B, ViridaA1, A3	P9325 M1900B+M 1196T	P8325 M1900B+ M1191T	P7325	P6325 M1900B+ M1192T	P5325 M1900B+ M1193T	P4325	P3325	3m	CN39	Direct 直连
MP20,MP30, MP40/50/60/70/90. VM4,VM6,VM8, M3046A	P9325A M1196A P9125A	P8325A M1191A/B M1191AL/BL P8125A	P7325A M1196A P7125A	P6325A M1192A P6125A	P5325A M1193A P5125A	P4325A	P3325A M1194A P3125A	3m 1m	CN40	Direct 直连 See cable
Philips-masim	P9325B	P8325B	P7325B	P6325B	P5325B	P4325B	P3325B	3m	CN41	Direct 直连
<b>Physio-Control</b>										
Lifepak 12	P9319	P8319	P7319	P6319	P5319	P4319	P3319	3m	CN30	Direct 直连
Lifepack 12, 20(Masim tech) 	P9315	P8315	P7315	P6315	P5315	P4315	P3315	3m	CN22	Direct 直连
<b>RGB</b>										
Omicrom FT, Omicron PREMIUM	P9303	P8303	P7303	P6303	P5303	P4303	P3303	3m	CP7	Direct 直连
<b>Schiller/ODM/ Brucker 席勒</b>										
CM8	P9326	P8326	P7326	P6326	P5326	P4326	P3326	3m	CN42	Direct 直连
TM7,Argus LCM	P9326A	P8326A	P7326A	P6326A	P5326A	P4326A	P3326A	3m	CN42	Direct 直连
TM 910	P9326B	P8326B	P7326B	P6326B	P5326B	P4326B	P3326B	3m	CN43	Direct 直连
Schiller Truscope 	P9326C	P8326C	P7326C	P6326C	P5326C	P4326C	P3326C	3m	CP7	Direct 直连
<b>Spacelabs-Burdick 太空</b>										
PRO2	P9303A	P8303A	P7303A	P6303A	P5303A	P4303A	P3303A	3m	CN04	Direct 直连
<b>Spacelabs 太空</b>										
90419, 90501, 90651A-03	P9323 369083-001	P8323	P7323	P6323	P5323	P4323	P3323	3m	CN36	Direct 直连
90351-0/6, 90465, 90466, 90467, 90651A-08, PC Express, 90489, IM77	P9327 015-0130-00+175-0646-00	P8327	P7327	P6327	P5327	P4327	P3327	3m	CN44	Direct 直连
90496 Ultraview	P9327A 015-0130-00 + 700-0029-00	P8327A	P7327A	P6327A	P5327A	P4327A	P3327A	3m	CN45	Direct 直连
(Masim tech)	P9327B	P8327B	P7327B	P6327B	P5327B	P4327B	P3327B	3m	CN48	Direct 直连
<b>Welch Allyn Protocol 伟伦</b>										
Propaq 206-EL, 204-EL, 202-EL, 106-EL	P9319	P8319	P7319	P6319	P5319	P4319	P3319	3m	CN30	Direct 直连
<b>Weinmann</b>										
Oxycount Mini	P9303A	P8303A	P7303A	P6303A	P5303A	P4303A	P3303A	3m	CN04	Direct 直连





Monitor Manufacturers 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能 Y 型	Ear clip 成人耳夹	Len. 长度	Conn. 接头	Adapter Cable 转接线
<b>Biolight 宝莱特</b>										
	P9305	P8305	P7305	P6305	P5305	P4305	P3305	3m	CP5	Direct 直连
Biolight M69; M9000 Biocare Neusoft NSC-M10; Contec CSM9000, CSM6000 (Digital)	P9305A	P8305A	P7305A	P6305A	P5305A	P4305A	P3305A	3m	CP5	Direct 直连
M700,M800 (Digital)	P9105C	P9105C	P9105C	P9105C	P9105C	P9105C	P9105C	1m	CN30	Direct 直连
M7000 (Digital)	P9305D	P8305D	P7305D	P6305D	P5305D	P4305D	P3305D	3m	CN61	Direct 直连
(Digital)	P9305E	P8305E	P7305E	P6305E	P5305E	P4305E	P3305E	3m	CN62	Direct 直连
(Nellcor Oximax)	P9305F	P8305F	P7305F	P6305F	P5305F	P4305F	P3305F	3m	CN61	Direct 直连
A8,M9000 (Digital)	P9305H	P8305H	P7305H	P6305H	P5305H	P4305H	P3305H	3m	CN82	Direct 直连
<b>Comen 科曼</b>										
CM8000A	P9361	P8361	P7361	P6361	P5361	P4361	P3361	3m	CP6	Direct 直连
(Digital)	P9361A	P8361A	P7361A	P6361A	P5361A	P4361A	P3361A	3m	CP6	Direct 直连
	P9361B	P8361B	P7361B	P6361B	P5361B	P4361B	P3361B	3m	CP6R	Direct 直连
(Digital)	P9361C	P8361C	P7361C	P6361C	P5361C	P4361C	P3361C	3m	CP6R	Direct 直连
Oximax	P9361D	P8361D	P7361D	P6361D	P5361D	P4361D	P3361D	3m	CP8R-R	Direct 直连
<b>Choice Electronic 超思</b>										
All(BCI tech)	P9303A	P8303A	P7303A	P6303A	P5303A	P4303A	P3303A	3m	CN81	Direct 直连
M3	P9342	P8342	P7342	P6342	P5342	P4342	P3342	3m	CP6	Direct 直连
(Digital)	P9342A	P8342A	P7342A	P6342A	P5342A	P4342A	P3342A	3m	CP6	Direct 直连
<b>Goldway 金科威</b>										
940	P9311	P8311	P7311	P6311	P5311	P4311	P3311	3m	CM5	Direct 直连
南京世帝;深迈 900F;科力威; 思创 7000C; 宏润达 EM6/EM7	P9311A	P8311A	P7311A	P6311A	P5311A	P4311A	P3311A	3m	CP5	Direct 直连
	P9311B	P8311B	P7311B	P6311B	P5311B	P4311B	P3311B	3m	CP7	Direct 直连
UT4000F (Digital)	P9311D	P8311D	P7311D	P6311D	P5311D	P4311D	P3311D	3m	CP5	Direct 直连
Nellcor oximax	P9311E	P8311E	P7311E	P6311E	P5311E	P4311E	P3311E	3m	CP7	Direct 直连
<b>Guoteng 珠海国腾</b>										
	P9350	P8350	P7350	P6350	P5350	P4350	P3350	3m	CP6	Direct 直连
<b>M&amp;B 北京麦邦</b>										
CD531, CD2000	P9351	P8351	P7351	P6351	P5351	P4351	P3351	3m	CN70	Direct 直连
<b>Neusoft 东软</b>										
	P9343	P8343	P7343	P6343	P5343	P4343	P3343	3m	CP7	Direct 直连
General 杰纳瑞 (Digital)	P9343A	P8343A	P7343A	P6343A	P5343A	P4343A	P3343A	3m	CN39	Direct 直连
<b>Solaris/New Tech 索莱瑞/纽泰克</b>										
(Digital)	P9320	P8320	P7320	P6320	P5320	P4320	P3320	3m	CP5	Direct 直连
General 杰纳瑞 (Digital)	P9320A	P8320A	P7320A	P6320A	P5320A	P4320A	P3320A	3m	CP6R	Direct 直连
(Digital)	P9120B	P8120B	P7120B	P6120B	P5120B	P4120B	P3120B	1m	CN30	See cable
世博达(Digital)	P9320C	P8320C	P7320C	P6320C	P5320C	P4320C	P3320C	3m	CP7	Direct 直连
科瑞康 PC-5000(Digital)	P9320D	P8320D	P7320D	P6320D	P5320D	P4320D	P3320D	3m	CP5T	Direct 直连

Monitor Manufacturers 监护仪厂家	Adult clip 成人指夹	Adult soft 成人软指套	Pedia. clip 儿童指夹	Pedia. Soft 儿童软指套	Neonate wrap 新生儿包裹	Multi-Y 多功能 Y 型	Ear clip 成人耳夹	Len. 长度	Conn. 接头	Adapter Cable 转接线
<b>Solaris/New Tech 索莱瑞/纽泰克</b>										
(Digital)	P9320E	P8320E	P7320E	P6320E	P5320E	P4320E	P3320E	3m	CN40	Direct 直连
	P9320F	P8320F	P7320F	P6320F	P5320F	P4320F	P3320F	3m	CP5	Direct 直连
<b>3F 上海三埃弗</b>										
	P9354	P8354	P7354	P6354	P5354	P4354	P3354	3m	CP6	Direct 直连
<b>Sunmind 深迈</b>										
	P9348	P8348	P7348	P6348	P5348	P4348	P3348	3m	CN40	Direct 直连
(Digital)	P9348A	P8348A	P7348A	P6348A	P5348A	P4348A	P3348A	3m	CN40	Direct 直连
<b>Zondan 和心重典</b>										
ZD120	P9345	P8345	P7345	P6345	P5345	P4345	P3345	3m	CP6	Direct 直连
<b>Jeteem 井田医电</b>										
BD9000	P9349	P8349	P7349	P6349	P5349	P4349	P3349	3m	CP6R	Direct 直连
<b>Adecon 爱德康</b>										
	P9378	P8378	P7378	P6378	P5378	P4378	P3378	3m	CP6	Direct 直连
<b>Perlong 南京普朗</b>										
	P9364	P8364	P7364	P6364	P5364	P4364	P3364	3m	CP6	Direct 直连

